# HAXLONE

## Cross-linked sodium hyaluronate viscosupplement 2.0%

**DESCRIPTION**

HAXL One is a biocompatible hydrogel of crosslinked sodium hyaluronate obtained by biofermentation and highly purified, dissolved in a physiological and isotonic buffer. The product is an injectable hydrogel, sterile and apyrogenic, with a low level of endotoxins.

**COMPOSITION**

HAXL One is formulated with highly purified sodium hyaluronate, pharmacopeia quality salts and water for injection (WFI):

COMPOSITION	(mg / mL)
Sodium hyaluronate cross-linked	20
Disodium phosphate dodecahydrate	0.6
Sodium dihydrogen phosphate dihydrate	0.05
Sodium chloride	8
Water For Injection (WFI)	q.s

**PACKAGING**

HAXL One is supplied in a graduated, pre-filled and disposable Type I glass sterile syringe with a Luer-lock system, containing 4.4 mL of solution.

**INTENDED USE**

HAXL One is intended to be injected intra-articularly in the synovial space of the knee as viscosupplementation medium for synovial fluid when its viscosity has been reduced due to osteoarthritis process in patients who have not responded adequately to conservative non-pharmacological therapy and simple analgesics.

**PROPERTIES**

HAXL One reduces the pain and increases the functionality of the affected joint. HAXL One provides a lubricant and shock-absorption effect on the affected joint. HAXL One improves the quality of life of the patient.

**ADMINISTRACION**

HAXL One must be administered by intra-articular injection by qualified medical practitioner following all the rules regarding the aseptic procedure and the injection technique. HAXL One should be injected with a suitable sterile needle (18-21G needles are recommended).

HAXL One is recommended to be administered in a single session. The amount of product to be administered is at the discretion of the specialist doctor. Discard any unused portion of the syringe.

**DIRECTIONS FOR USE**

- Prior to the injection, the area to be treated should be disinfected thoroughly and bring the product to room temperature.
- Peel off the sterile barrier and remove the syringe from inside.
- Remove the tip cap of the syringe without touching the syringe tip, in order to maintain it sterile.
- Attach a suitable (18-21G recommended) sterile needle to the tip of the syringe. Verify that the needle is properly locked to theLuer -lock adapter of the syringe before using it.
- Remove air from the syringe, if necessary.
- Inject intra-articularly into the synovial space of affected knee, following aseptic injection procedures.
- After finishing the treatment session, discard the syringe and the remaining product. The amount of product to be administered is at the discretion of the specialist doctor; the recommended dose is 3mL to 4mL for each knee joint.

**CONTRAINDICATIONS**

Do not administer HAXL One to patients with known hypersensitivity (allergy) to hyaluronic acid and / or other components of the product. Do not administer HAXL One to patients with acute or chronic dermatological disorders or an inflammatory and / or infectious type skin problem.

Do not administer HAXL One to patients with coagulation defects and / or who are treated with anticoagulant medication.

HAXL One should not be administered to pregnant or breastfeeding women. HAXL One should not be administered to children.

The use of HAXL One in any joint other than the knee and for other conditions has not been studied.

**WARNINGS AND PRECAUTIONS**

HAXL One must be injected strictly into the intra-articular space following aseptic injection procedures.

HAXL One should be administered with sterile Luer-lock needles.

Do not inject into a blood vessel or surrounding tissues. To avoid pain after the injection, a recommendation should be given to keep the local site still after the injection.

The product should not be used if the package is damaged or broken. The product should not be used after the expiration date.

The product is for single use. No reuse; risk of infection.

Do not re-sterilize, re-sterilization can affect the physico-chemical properties of the product and compromise its effectiveness and safety.

The syringe and the used needle should be discarded in a container provided for this purpose.

Keep out of children’s reach.

**INTERACTIONS**

Do not use with other products containing quaternary ammonium salts as benzalkoni-um chloride or with any surgical material treated with this type of substances because sodium hyaluronate can precipitate in their presence. In order to prevent any possible interaction, avoid joint administration with other intra-articular products.

There is no known interaction with other local anesthetics.

**ADVERSE EVENTS**

Viscoelastic supplementation with hyaluronic acid in joints is safe. Moderate episodes of arthralgia, joint swelling or stiffness, and transient discomfort and / or pain at the injection site following intra-articular injection have occasionally been observed. To a lesser extent, there may be some case of joint effusion or minor local skin conditions such as erythema or rash.

**STORAGE**

Store between 5°C and 25°C. Protect from light and avoid freezing.

**SELL BY DATE**

Do not use after the date marked on the package.

**MANUFACTURED BY & LEGAL MANUFACTURER**

	i+Med S.Coop Parque Tecnológico de Álava Calle de Albert Einstein, 15. P. 15 01510 Vitoria-Gasteiz (Álava - Spain) <b>www.imasmed.com</b>
---	---

Segregation and waste management according to current legislation.

												
---	---	---	---	---	--	---	---	---	---	---	---	---

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