
EN
NL
FR
DE
ITA
PL

PRP
TUBES

This Instructions for Use document covers all brands of the product with reference 1003&1005.

The name or trade name of the medical device

PRP TUBES

The product is a sterile, single use medical device.

The name and the address of the manufacturer

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PRP TUBES

Intended use of the product

The device is designed to be used for the safe and rapid preparation of autologous platelet-rich-plasma (PRP) from a small sample of blood at the patient's point of care. The PRP is prepared prior to application to a local area of soft tissue or hard tissue (including bone) site as deemed necessary by the clinical use requirements. And its intended use is for the application of PRP in order to achieve the stimulation and acceleration of soft tissue and hard tissue (including bone) healing.

Intended patient population

Platelet rich plasma is an autologous application technique therefore patient applications are inherently safe from infectious diseases such as HIV, Hepatitis, Creutzfeldt-Jakob disease.

This application can be used safely and effectively on all patients but the safety and effectiveness of PRP treatment have not been established in children and pregnant women.

The device is intended to be used for the treatment of the following conditions:

- alopecia areata,
- joint pain (particularly in patients with osteoarthritis of the knee),
- rotator cuff injuries,
- and lateral epicondylitis (tennis elbow).

The performance of the device intended by the manufacturer

PRP (Platelet Rich Plasma) is used for many years in the frame of regenerative medicine, which allows the cellular healing on bone and soft tissue recovery. Bone and soft tissue healing occurs by the re-injection of the platelet concentration locally in the related part of the body.

Currently accepted technique for avoiding coagulation of the whole blood, referred to is Sodium Citrate. The procedure uses anticoagulated blood, collected by routine phlebotomy. Buffered Citrate Solution blocks the coagulation by mixing gently the whole blood together with 3,20% concentrated 0,1M of Sodium Citrate. The prepared material is then centrifuged to

obtain PRP. The recommended G Force for concentration is given in the "Procedure" section according to the centrifuge rotor.

The blood separation system makes use of the advantage of the effect of centrifugation force on the weight of the cells. The separation occurs during centrifugation when the denser blood components move due to a gravity force.

PRP TUBES should not be shaken after the centrifugation. Depending on the centrifuge rotor used and the desired protocol, PRP is harvested carefully above the erythrocytes (Red Blood Cell layer).

The result is a convenient system, suitable for collection of the whole blood and concentration of Platelet Rich Plasma. PRP Tubes are interior sterilized using irradiation. PRP TUBES reduces the risk of air contamination, and eliminates the need for any other tubes or pipettes.

Clinical safety information

The safety and efficacy of this device has been proven by Risk Analysis Report and Clinical Evaluation Report, which are attached to the Technical File of the PRP TUBES.

Residual risk and side effects information

According to the risk assessment, risk management, risk report procedures of this device, it shows no residual risks. With respect to the GSPR 3, the risk management report showed that there were no serious side effects and the performance of the product was confirmed by the clinical evaluation report.

Intended user

This product should be used by physicians or physician directed allied health care professionals with adequate training in the use of the device.

Medical devices to be used with the device

- Blood Collection Needle / Phlebotomy Set / Butterfly Set
 - Phlebotomy set / butterfly set / blood collection needle shall not be thinner than 21G.
- Centrifuge
 - The centrifuge should be set as specified in the "Procedure" section.
 - PRP TUBES have dimensions of 16X100 mm. Centrifuge bucket shall be selected accordingly.
 - Place the tubes in the centrifuge so that they are of equal weight.
 - Avoid using unbalanced centrifugal devices.
- Long Needle and Syringe
 - Long needle shall be 88mm-90mm 18G-19G.
 - Long needle used to collect PRP layer from the tube is not intended to be used for the injection procedure.
 - Syringe used to collect the PRP layer shall have a minimum capacity of 4ml and be a luerlocked model.

Procedure

Apply the following instructions to obtain Platelet Rich Plasma by using PRP TUBES:

- Open the box.
- Open the blister package.
- Apply the blood collection device to sanitized part of the patient's arm. (Venous)
- PRP TUBES is a vacuumed tube. The tube collects enough amount of blood and stops

itself.

- Centrifuge the PRP Tubes in a centrifuge according to the procedures outlined below (The doctor must choose the procedure).

	Fixed Rotor Centrifuge	Swing Rotor Centrifuge
Pure PRP Protocol	500 RCF & 3 Minute	1150 RCF & 4 Minute
Amber PRP Protocol	500 RCF & 6 Minute	1150 RCF & 8 Minute
500 RCF=2000 RPM, 1150 RCF=2900 RPM		
AMBER PRP: LR-PRP (Leukocyte Rich PRP) PURE PRP: LP-PRP (Leukocyte Poor PRP)		

- Buckets shall be compatible for PRP Tubes (16 X 100 mm tube size).
- Attach a long needle to 5ml Luer-lok Syringe.
- Collect the PRP by sticking the Long needle into the PRP Tube above the RBC Layer (Do not open the Lid of the tube unless it is performed under Laminar Flow Conditions)
- After the Pure PRP Protocol, collect all the plasma above the erythrocyte layer
- After the Amber PRP Protocol, collect the buffy coat along with 2ml plasma above the erythrocyte layer.
- Carefully collect the PRP according to the selected protocol
- Apply the injection method to the patient.

Storage, shelf life and handling information

- Do not expose to direct sunlight.
- Do not use if the inner blister package is damaged or opened.
- Expiration date is indicated on the label.
- Storage temperature is between 5°C and 30°C

Precautions, warnings, measures to be taken and limitations

- Use proper safety precautions to guard against needle sticks.
- Follow manufacturer's instructions when using centrifuge.
- Do not use sterile component of this kit if package is opened or damaged.
- Do not open the lids of the tubes.
- Single use device, do not reuse.
- The user should be thoroughly familiar with the equipment and the surgical procedure prior to using this device.
- The patient should be made aware of general risks associated with treatment and possible adverse effects.
- Do not dispose of this product in the unsorted municipal waste stream.
- Long Needle is not intended to be used on humans. It is only intended to be used to collect the PRP Layer inside from the PRP Tube(s).

Possible adverse events

- Damage to blood vessels, hematoma / delayed wound healing and/or infection associated with blood draw.
- Temporary or permanent nerve damage that may result in pain or numbness.
- Early or late postoperative infection.

Possible contraindicative situations

Platelet Dysfunction Syndrome, Critical Thrombocytopaenia, Hypofibrinogenaemia, Hemodynamic instability, Auto-immune diseases, Malignancy, Sepsis, Acute Chronic infections, Chronic Liver Pathology, Anti-coagulation Therapy (ASA therapy etc.), Pregnancy, Newborns, Not recommended for use under 18 years of age.

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